

# LACTOTEST 202

## Hydrogen, methane and carbon dioxide breath test monitor


The **Lactotest 202** has been developed to provide quick and reliable diagnosis of gastro-intestinal disorders measuring carbohydrate breakdown deficiencies and/or malabsorption. Particularly efficient for lactose intolerance screening, the **Lactotest 202** also gives good indication of small intestinal bacterial overgrowth (SIBO) and evaluates intestinal transit time. The **Lactotest 202** range is composed of:

### Lactotest 202 "Xtend" (H<sub>2</sub> + CO<sub>2</sub> + CH<sub>4</sub>) & Lactotest 202 "Control" (H<sub>2</sub> + CO<sub>2</sub>)

- The electrochemical sensor with thermal compensation for H<sub>2</sub> and infra-red analyzers for CH<sub>4</sub> and CO<sub>2</sub> deliver an accurate test measurement in less than one minute. The Lactotest 202 Control can be upgraded to CH<sub>4</sub> detection.



Lactotest 202- Xtend



- Approximately 35% of the population produces methane, alternatively or in addition to hydrogen\*. The role of CH<sub>4</sub> is to ensure that these patients are not miss-diagnosed, therefore reducing the percentage of false negative tests.
- The CO<sub>2</sub> concentration is virtually constant in the alveolar sample. If the sample is not alveolar, the results are underestimated. With the CO<sub>2</sub> correction factor (applied to H<sub>2</sub> and CH<sub>4</sub>), the **Lactotest 202** implements an essential verification of the results' reliability.
- The technique of exhaled breath monitoring is usually well tolerated by patients of all ages. The test is easy, **non-invasive** and performed after a short period of fasting.
- The **Lactotest 202** works with an automatic pump for test sampling and flushing of residual gas, making it immediately ready for the next test.

\* Clinical Gastroenterology Hepatology, 2006 Feb; 4 (2): 123-190



## Patient/User benefits

<b>Applications</b>	Carbohydrate malabsorption detection Lactose malabsorption and/or intolerance Fructose/Glucose/Lactulose/Sorbitol/Sucrose/d-Xylose/Xylitol intolerance Small Intestinal Bacterial Overgrowth (SIBO) diagnosis Transit time diagnosis
<b>Non invasive</b>	Sampling by exhaled air collection
<b>Software</b>	Patient Data Interface (PDI): patients' database integration; import/export of results into Hospital Information Systems (HIS); report creation (MS Word, pdf, etc.)
<b>Calibration</b>	Once a day with pre-mixed calibration gas

## How to choose the Lactotest 202 you need?

	LACTOTEST 202 « CONTROL »	LACTOTEST 202 « XTEND »
ELECTROCHEMICAL H2 ANALYSER		
INFRA RED CH4 ANALYSER		
INFRA RED CO2 ANALYSER		
INJECTION / FLUSHING / SAMPLING	MANUAL / AUTOMATIC / SYRINGE	AUTOMATIC / AUTOMATIC / BAG

## Specifications

<b>Measurement range</b>	H2: 0-200 ppm; CH4: 0-500 ppm; CO2 0-10%
<b>Accuracy</b>	± 2 % fs
<b>Detection principle</b>	H2: Electrochemical sensor with thermal compensation CH4 & CO2: Infra-red analysers
<b>Sensor operating life (and warranty)</b>	H2: minimum 3 years (1 year warranty) CH4 & CO2: 10 years (1 year warranty)
<b>Sensor resolution</b>	H2: 1 ppm CH4: 1 ppm & CO2: 0.01%
<b>Warm-up time</b>	202 Control: 30 min.; 202 Xtend: 60 min.

Certification / Safety standards  
93/42/EEC & 2007/47/EC Medical Devices Directive

## MEC, YOUR PARTNER FOR THE LONG TERM ...



Medical Electronic Construction R&D sprl  
Rue Prévinair 64 Previnairstraat  
Bruxelles 1070 Brussel  
BELGIUM

**Phone** + 32(0) 2 558.00.60  
**Fax** + 32(0) 2 558.00.69  
**e-mail** info@mecrd.eu  
**Web** www.mecrd.eu

